

SPECTROFLUOROMETER

Brand-Model

FluoroMax-4, HORIBA

Descriptions

Short description of fungsi dan spec (apa yang special?)

Service Offer??

Type of samples??

Water Raman signal-to-noise ratio: 10,000:1 (FSD method)

Photon Counting for Ultimate Sensitivity

Fast Scanning capability - up to 80 nm/second

Specifications

Slit: Continuously variable from 0 to 30nm

Minimum step: 0.0525nm

Monochromators: All reflective optics, Czerny-Turner spectrometers

Accuracy: 0.5nm

Service location

03-01-10 *Material Characterization Laboratory*

Block: T05

Further Information

- Assoc. Prof. Dr. Abd Rahman Tamuri
- E-mail : rahmantamuri@utm.my

Operator

- Ms Rahizah Binti Abd Rahman
- E-mail : rahizah.ar@utm.my

Notes:

